



Bainmount Pneumatic

Pneumatic Specimen Mounting Press



Mounting the specimen after sectioning is often necessary for subsequent handling and Metallographic polishing.

Bainmount Pneumatic



Mounting of the specimen after sectioning is often necessary for subsequent handling and Metallographic polishing. Mounting has several benefits, especially in hand polishing when specimen flatness and edge retention are important.

Specifications

Corrosion resistant steel Cabinet

Pneumatic operation System

Digital Temperature Indicator and digital timer with LCD Display

Quick Ram return to speed up reloading

1600 Watts mould Heater

Dovetail slot for easy changeability of Mould Assemblies

Water Cooling System (manual)

Timer-Buzzer indication for heating and cooling cycles

Wide range of sample sizes can be accommodated

Provision to accommodate 1 ½" & 2" Mould Sizes

1¼" comes with the machine as a standard

Power Supply : 230V / 50Hz (1 + N + E)

Machine Size: 580x620 x 590mm

Optional Accessories

- Mould Assembly 1½"
- Mould Assembly 2"
- Recirculation Water Cooling System 50 Ltrs capacity
- 1 ¼" Twin Mould design
- 1 ½" Twin Mould design
- 2" Twin Mould design
- Silent Air compressor

MPP006

Chennai Metco Pvt Ltd

SP100A, Ambattur Industrial Estate,
Chennai 600 058, India.

Ph: 044-2688 0005, 044-3999 5000, Fax: 044-2688 0850

e-mail: mktg@chennaietco.com, mktg@chennaietco.co.in

www.chennaietco.com

