



TEMPERATURE: 157.0  
HEATING TIME: 00.01  
AMPS: 2.0A

# Bainmount P Auto

Pneumatic - Automatic Specimen  
Mounting Press

Digital Temperature  
indicator and Digital  
Timer with Display.  
Automatic water  
cooling system.



[www.chennaietco.com](http://www.chennaietco.com)

# Bainmount P Auto



Mounting of the specimen after sectioning is often necessary for subsequent handling and Metallographic polishing. Mounting has several benefits, especially in hand polishing when specimen flatness and edge retention are important.

## Specifications

- Steel Cabinet with built in Pneumatic System
- Digital Temperature indicator and Digital timer with LCD Display
- Quick Ram return to speed up reloading and automatic cycle system
- Automatic Water Cooling System
- Thermostatically controlled 1600 Watts Mould heater
- Buzzer indication for cooling cycles
- Wide Range of Sample size can be accommodated
- With Interchangeable 1¼" , 1½" 2" mould sizes
- 1¼" comes with machine as a standard
- Power supply : 230V /50Hz ( 1 + N + E )
- Machine Size:580 x 630x 600 mm
- Machine Color: Enclosure RAL 9003, Base RAL 7021

## Optional Accessories

- Mould Assembly 1½"
- Mould Assembly 2"
- Recirculation Water Cooling System 50 Ltrs capacity
- 1 ¼" Twin Mould design
- 1 ½" Twin Mould design
- 2" Twin Mould design
- Silent Air compressor

MPP006A

## Chennai Metco Pvt Ltd

SP100A, Ambattur Industrial Estate,  
Chennai 600 058, India.

Ph: 044-2688 0005, 044-3999 5000, Fax: 044-2688 0850

e-mail: mktg@chennaietco.com, mktg@chennaietco.co.in

[www.chennaietco.com](http://www.chennaietco.com)

